

NOSOFLOOR®

Highly concentrated disinfectant
for floors and surfaces

NOSOFLOOR is a high efficacy concentrated solution for the cleaning and disinfection of hard surfaces as walls, floors, toilet bowls, door handles, etc. It combines a broad spectrum of microbicidal properties and a very good cleaning power.

It is suitable as a routine disinfection procedure in operating rooms, intensive care units, corridors and patient-near areas.

NOSOFLOOR has an excellent compatibility with most materials.

It does not contain phenols, aldehydes, chlorine or EDTA.

Properties

- Bactericidal, fungicidal, tuberculocidal, mycobactericidal
- Virucidal (HBV, HIV, HCV, Herpes, Vaccinia, BVDV, Influenza, Ebola, Coronavirus...)
- Excellent cleaning properties
- Does not affect the medical equipment
- Leaves no residue after drying
- Does not contain phenols, aldehydes, chlorine or EDTA

Disinfecting properties

ACTIVITY SPECTRUM	STANDARD	STRAINS	CONTACT TIME
BACTERICIDAL (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	0.25% - 5 Min.
FUNGICIDAL (Dirty conditions)	EN 13624	Candida Albicans	0.25% - 5 Min.
VIRUCIDAL (Dirty conditions)	DVV ⁽¹⁾ /RKI ⁽²⁾ 2014	BVDV, Vaccinia, HBV, HIV, HCV, Ebola, Herpes, Influenza H1N1, H5N1, Coronavirus	1% - 15 Min. 0.75% - 30 Min
TUBERCULOCIDAL (Dirty conditions)	EN 14348	Mycobacterium Terrae (Surrogate M. tuberculosis)	2% - 15 Min. 1% - 60 Min.
MYCOBACTERICIDAL (Dirty conditions)	EN 14348	Mycobacterium Terrae Mycobacterium Avium	1% - 60 Min.

(1) DVV: Deutsche Vereinigung zur Bekämpfung der Viruskrankheiten (German Association for the Control of Virus Diseases)

(2) RKI: Robert Koch Institute - German Federal Health Authority

Compatibility

NOSOPROTECT is compatible with most materials such as stainless steel, aluminum, glass, ceramics, hard plastics, linoleum, ebonite, etc ...

Certifications

- Biocidal product according to Directive 98/8/CE

NOSOFLOOR is manufactured in the EU.

MEDALKAN satisfies the requirements of ISO 9001:2015 for quality management system and the requirements of ISO 13485:2016 for the design and manufacture of medical devices.

For further informations, please visit www.medalkan.gr

MEDALKAN - 102, Michalakopoulou street, 11528 Athens - Greece - Tel.: (+30) 210-7484847 - Fax: (+30) 210-7772009 - www.medalkan.com - contact@medalkan.gr



Packaging

- 5 litres canister. (Ref. 20025)

Physical properties

- Appearance: Transparent blue solution
- Density: 0.99 g/cm³ at 20°C
- pH: 12.0-12.8 at 20°C
- pH (1%): 9.5-10.5 at 20°C
- Odour: Neutral
- Storage: 5°C - 35°C
- Stability: 3 Years
- Biodegradability: According to OCDE 301D

Composition

N-(3-aminopropyl)-N-dodecylpropano-1,3-diamine, didecyl-dimethyl ammonium chloride, non-ionic surfactants <5%, isopropyl alcohol, corrosion inhibitor, excipients

Instructions for use

Dilute NOSOFLOOR in cold water according to the desired disinfecting activity (See table). Stir gently. Make sure to wet surfaces completely and keep them wet for the whole exposure time. A rinsing procedure is required only for surfaces in contact with the skin.

Do not mix with other disinfectants.

The packaging must be closed tightly after use.

Store the product between + 5 ° C and + 35 ° C.

Use biocides with care. Before use, read the label and product information.



Medalkan
HYGIENE & DISINFECTION