

MEDASEPT® 100

Fast acting broad spectrum surface disinfectant



MEDASEPT 100 is a fast acting and broad spectrum sprayable disinfectant for the medical devices' surfaces after each session.

It is specifically designed to avoid cross contamination between patients. It is recommended for the disinfection of medical equipment surfaces (benches, unit, handles, spittoon, etc...) that are in direct contact with patients and medical staff.

It has excellent disinfecting properties and leaves no traces after drying. MEDASEPT 100 has a mild odor and does not affect the medical equipment.

Properties

- Bactericidal, fungicidal, tuberculocidal, mycobactericidal
- Virucidal (HBV, HIV, HCV, Herpes, Vaccinia, BVDV, Influenza, Ebola, Coronavirus...)
- Active in 30 seconds
- Does not affect the medical equipment
- Leaves no residue after drying
- Does not contain phenols, aldehydes, chlorine or EDTA

Disinfecting properties

ACTIVITY SPECTRUM	STANDARD	STRAINS	CONTACT TIME
BACTERICIDAL (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	30 sec.
FUNGICIDAL (Dirty conditions)	EN 13624	Candida Albicans Aspergillus Brasiliensis (Fungicidal)	30 sec. 5 min.
VIRUCIDAL (Dirty conditions)	DVV ⁽¹⁾ /RKI ⁽²⁾ 2014	BVDV, Vaccinia, HBV, HIV, HCV, Ebola, Herpes, Influenza H1N1, H5N1, Coronavirus	30 sec.
TUBERCULOCIDAL (Dirty conditions)	EN 14348	Mycobacterium Terraе (Surrogate M. tuberculosis)	3 min.
MYCOBACTERICIDAL (Dirty conditions)	EN 14348	Mycobacterium Terraе Mycobacterium Avium	3 min.

(1) DVV: Deutsche Vereinigung zur Bekämpfung der Viruskrankheiten (German Association for the Control of Virus Diseases)

(2) RKI: Robert Koch Institute - German Federal Health Authority

Certifications

- CE mark according to the medical devices Directive (Directive 93/42/EEC)
- Medical device class IIa

MEDASEPT 100 is manufactured in Greece.

MEDALKAN satisfies the requirements of ISO 9001 for quality management system and the requirements of ISO 13485 for the design and manufacture of medical devices.

For further informations, please visit www.medalkan.com

MEDALKAN - Michalakopoulou 102, 115-28 Athens - Tel: (+30) 210-7484847 - Fax: (+30) 210-7772009 - www.medalkan.com - contact@medalkan.gr



Packaging

- One litre bottle with spray (Ref. 20001)
- 5 litres canister (Ref. 20002)

Physical properties

- Appearance: Transparent solution
- Density: 0.97 g/cm³ at 20 °C
- pH: 9.0-9.6 at 20 °C
- Odour: Mild (alcohol)
- Storage: 5 °C - 35 °C
- Stability: 3 Years
- Biodegradability: According to OCDE 301D

Instructions for use

Spray MEDASEPT 100 directly on the surface to be cleaned. Let the solution act for at least 30 seconds and wipe with a single-use paper. Let dry. Do not rinse. Do not mix with other disinfectants. Store between + 5 ° C and + 35 ° C.

Composition

Isopropyl alcohol, didecyl-dimethyl ammonium chloride, N- (3-aminopropyl) -N dodecylpropano-1,3-diamine, excipients.



Medalkan
HYGIENE & DISINFECTION